

PRODUCTS CATALOGUE 2017

ORTHOPAEDIC PRODUCTS

REHABILITATION PRODUCTS

ASSISTIVE DEVICES





ARmedical®

CONSUMER LIDER OF QUALITY 2016

We are pleased to announce that ARmedical brand received award
CONSUMER LEADERSHIP OF QUALITY 2016




Orthopaedic pillows

MFP-5030**Profiled orthopaedic pillow
"EXCLUSIVE DREAM"**

Made from an innovative thermoforming foam with a shape "memory" effect. This allows the pillow to memorise a shape of the head and as a result lower the pressure on the body keeping proper blood circulation.

- » Provides optimal alignment of the cervical vertebrae during sleep.
- » Thanks to the ergonomic profiling, pillow corrects the position of the spine during sleep through the support of the head, neck and shoulders.
- » Relieves spinal, neck and shoulder pain by reducing pressure points.
- » Innovative foam prevents the growth of microorganisms, allergens and absorption of dust particles.

**Highly recommended in the case of:**

- » spine pain
- » shoulder pain
- » hand numbness
- » snoring
- » insomnia
- » nervous tensions

Pillow has a super soft, innovative micro-fibre pillowslip (zippered) which allows washing of the pillow at the same time preserving its aesthetic appearance.

Dimensions: 50 x 30 x 10 cm; **Cover:** micro-fiber, **Colour:** Ecrú
3 years WARRANTY

Cover available separately: cat. C-5030

ARmedical®
Pillow Line



MFP-5030MV

Profiled orthopaedic pillow with massage & ventilation "EXCLUSIVE DREAM"

Made from an innovative thermoforming foam with a shape "memory" effect. This allows the pillow to memorise a shape of the head and as a result lower the pressure on the body keeping proper blood circulation.

MFP-5030 orthopedic pillow is used to prevent and ease the disorders and to mitigate or compensate the effects of cervical spine injury by ensuring its proper alignment during sleep.

Highly recommended in the case of:

- » spine pain
- » shoulder pain
- » hand numbness
- » snoring
- » insomnia
- » nervous tensions

Properties:

- » Relieves spinal, neck and shoulder pain by reducing pressure points
- » Guarantees correct positioning of the cervical spine during sleep through the proper support of the head, neck and shoulders
- » Prevents morning stiffening "numbness" of hands
- » Perfectly adapts to the anatomical shape of the body
- » Prevents the growth of microorganisms, allergens and absorption of dust particles

Pillow has a super soft, innovative micro-fiber pillowslip (zippered) which allows washing of the pillow at the same time preserving its aesthetic appearance.

Dimensions: 50x30x10cm; **Cover:** micro-fiber, **Colour:** Ecu; 3 years WARRANTY

Cover available separately: cat. C-5030



ZARmedical
Pillow Line

**MFP-5535****Profiled orthopaedic pillow "EXCLUSIVE DREAM"**

Made from an innovative thermoforming foam with a shape "memory" effect. This allows the pillow to memorise a shape of the head and as a result lower the pressure on the body keeping proper blood circulation.

Anatomically shaped cushion adapts perfectly to the shape of the patient's body whether the patient is sleeping on his/her side or on the back. The head almost submerges in the pillow's foam, thereby increasing the contact surface of the body/head with the surface as a result reducing pressure force.

Highly recommended in the case of:

- » spine pain
- » shoulder pain
- » hand numbness
- » snoring
- » insomnia
- » nervous tensions

Pillow has a super soft, innovative micro-fibre pillowslip (zippered) which allows washing of the pillow at the same time preserving its aesthetic appearance.

Dimensions: 55x35x12cm; **Cover:** micro-fibre; **Colour:** Ecru

3 years WARRANTY

Cover available separately: cat. C-5535

**MFP-5030BF****Profiled orthopaedic pillow "BAMBOO DREAM"**

Made from an innovative thermoforming foam with a shape "memory" effect. Innovative foam with the addition of active bamboo charcoal stimulates blood circulation, reduces stress and muscle tension. Active bamboo charcoal in the foam as well as natural bamboo fibre in a pouch prevent the growth of microorganisms, allergens and increase absorption of dust particles.

Thanks to the ergonomic profiling, pillow corrects the position of the spine during sleep through the support of the head, neck and shoulders. Relieves spinal, neck and shoulder pain by reducing pressure points. With thermoforming foam it perfectly adapts to the anatomical shape of the body.

Highly recommended in the case of:

- » spine pain
- » shoulder pain
- » hand numbness
- » snoring
- » insomnia
- » nervous tensions

Pillow comes with two covers, inner and outer. Inner pouch provides additional protection of the foam. Outer cover is made from innovative material with the addition of natural bamboo fibre. Provides proper air circulation without causing excessive sweating of the body. Comes with a zipper which allows washing of the pillow at the same time preserving its aesthetic appearance.

Properties:

- » Provides optimal alignment of the cervical vertebrae during sleep
- » Thanks to the ergonomic profiling, pillow corrects the position of the spine during sleep through the support of the head, neck and shoulders
- » Relieves spinal, neck and shoulder pain by reducing pressure points
- » Innovative foam with the addition of active bamboo charcoal prevents the growth of microorganisms, allergens and absorption of dust particles

Dimensions: 50x30x10cm; **Cover:** Jersey with natural bamboo fibre; **Colour:** Ecru

3 years WARRANTY

Cover available separately: cat. C-5030BF





MFP-5030BV Profiled orthopaedic pillow with massage & ventilation "Bamboo Dream"

Made from an innovative thermoforming foam with a shape "memory" effect. Innovative foam with the addition of active bamboo charcoal stimulates blood circulation, reduces stress and muscle tension. Active bamboo charcoal in the foam as well as natural bamboo fibre in a cover prevent the growth of microorganisms, allergens and increase absorption of dust particles.

- » Provides optimal alignment of the cervical vertebrae during sleep.
- » Thanks to the ergonomic profiling, pillow corrects the position of the spine during sleep through the support of the head, neck and shoulders.
- » Perfectly adapts to the anatomical shape of the body.
- » Relieves spinal, neck and shoulder pain by reducing pressure points.
- » Innovative foam with the addition of active bamboo charcoal prevents the growth of microorganisms, allergens and absorption of dust particles.

Ventilation holes provide proper air circulation without causing excessive sweating of the body increasing comfort at the same time. Due to its ergonomic profiling, pillow corrects the position of the spine during sleep through the support of the head, neck and shoulders. Relieves spinal, neck and shoulder pain by reducing pressure points. With thermoforming foam it perfectly adapts to the anatomical shape of the body.

Highly recommended in the case of:

- » spine pain
- » shoulder pain
- » hand numbness
- » snoring
- » insomnia
- » nervous tensions



Pillow comes with two covers, inner and outer. Inner cover provides additional protection of the foam. Outer cover is made from innovative material with the addition of natural bamboo fibre. Provides proper air circulation without causing excessive sweating of the body. Comes with a zipper which allows washing of the pillow at the same time preserving its aesthetic appearance.

Dimensions: 50x30x10cm; **Cover:** Jersey with natural bamboo fibre; **Colour:** Ecrú

3 years WARRANTY

Cover available separately: cat. C-5030BF



MFP-5232 Profiled orthopaedic pillow "BAMBOO DREAM"

Made from an innovative thermoforming foam with a shape "memory" effect. Innovative foam with the addition of active bamboo charcoal stimulates blood circulation, reduces stress and muscle tension. Active bamboo charcoal in the foam as well as natural bamboo fibre in a cover prevent the growth of microorganisms, allergens and increase absorption of dust particles.

Thanks to the ergonomic profiling, pillow corrects the position of the spine during sleep through the support of the head, neck and shoulders. Relieves spinal, neck and shoulder pain by reducing pressure points. With thermoforming foam it perfectly adapts to the anatomical shape of the body.

- » Provides optimal alignment of the cervical vertebrae during sleep.
- » Thanks to the ergonomic profiling, pillow corrects the position of the spine during sleep through the support of the head, neck and shoulders.
- » Perfectly adapts to the anatomical shape of the body.
- » Relieves spinal, neck and shoulder pain by reducing pressure points.
- » Innovative foam with the addition of active bamboo charcoal prevents the growth of microorganisms, allergens and absorption of dust particles.



Highly recommended in the case of:

- » spine pain
- » shoulder pain
- » hand numbness
- » snoring
- » insomnia
- » nervous tensions



Pillow comes with two covers, inner and outer. Inner cover provides additional protection of the foam. Outer cover is made from innovative material with the addition of natural bamboo fibre. Provides proper air circulation without causing excessive sweating of the body. Comes with a zipper which allows washing of the pillow at the same time preserving its aesthetic appearance.

Dimensions: 52x32x12cm; **Cover:** Jersey with natural bamboo fibre; **Colour:** Ecrú

3 years WARRANTY

Cover available separately: cat. C-5232



MFP-5030AF Profiled orthopaedic pillow "AloeVera Dream"

Made from an innovative thermoforming foam with a shape "memory" effect. Innovative foam with the addition of natural aloe vera extract stimulates blood circulation, reduces stress and muscle tension. Natural aloe vera extract in the foam as well as in a cover prevents the development of microorganisms, fungi and allergens.

Due to its natural aloe vera extract the pillow supports body regeneration during sleep and increases the comfort.

Ergonomic profiling corrects the position of the spine during sleep through the support of the head, neck and shoulders. Relieves spinal, neck and shoulder pain by reducing pressure points.



Highly recommended in the case of:

- » spine pain
- » shoulder pain
- » hand numbness
- » snoring
- » insomnia
- » nervous tensions

Dimensions: 50x30x11cm;
Cover: Jersey with natural aloe vera extract; **Colour:** Ecru

3 years **WARRANTY**

Cover available separately:
 cat. C-5030AF

Properties:

- » Due to its beneficial effect on health, aloe vera has been valued for many years in natural medicine
- » It naturally supports the regeneration of the body during sleep and increases its comfort
- » Due to its ergonomic profiling, pillow corrects the position of the spine during sleep through the support of the head, neck and shoulders
- » Relieves spinal, neck and shoulder pain by reducing pressure points
- » With thermoforming foam it perfectly adapts to the anatomical shape of the body
- » Innovative foam with the addition of natural aloe vera extract stimulates blood circulation, reduces stress and muscle tension
- » Prevents the development of microorganisms, fungi and allergens
- » Innovative, soft, breathable pillowslip with natural aloe vera extract comes with a zipper thus maintaining the pillow in proper hygienic condition is not a problem



Properties:

- » Due to its beneficial effect on health, aloe vera has been valued for many years in natural medicine
- » It naturally supports the regeneration of the body during sleep and increases its comfort
- » Due to its ergonomic profiling, pillow corrects the position of the spine during sleep through the support of the head, neck and shoulders
- » Relieves spinal, neck and shoulder pain by reducing pressure points
- » With thermoforming foam it perfectly adapts to the anatomical shape of the body
- » Innovative foam with the addition of natural aloe vera extract stimulates blood circulation, reduces stress and muscle tension
- » Prevents the development of microorganisms, fungi and allergens
- » Innovative, soft, breathable pillowslip with natural aloe vera extract comes with a zipper thus maintaining the pillow in proper hygienic condition is not a problem



MFP-5535AF Profiled orthopaedic pillow "AloeVera Dream"

Made from an innovative thermoforming foam with a shape "memory" effect. Innovative foam with the addition of natural aloe vera extract stimulates blood circulation, reduces stress and muscle tension. Natural aloe vera extract in the foam as well as in a cover prevent the growth of microorganisms, fungi and allergens.

Due to its natural aloe vera extract the pillow supports body regeneration during sleep and increases the comfort.

Ergonomic profiling corrects the position of the spine during sleep through the support of the head, neck and shoulders. Relieves spinal, neck and shoulder pain by reducing pressure points.



Highly recommended in the case of:

- » spine pain
- » shoulder pain
- » hand numbness
- » snoring
- » insomnia
- » nervous tensions

Dimensions: 55x35x12cm;
Cover: Jersey with Natural aloe vera extract;
Colour: Ecru

3 years **WARRANTY**

Cover available separately:
 cat. C-5535AF





MFP-5433EP

Double-sided profiled orthopaedic pillow "Ergonomic Dream"

ergonomic

Made from an innovative thermoforming foam with a shape "memory" effect. As a result, it perfectly adapts to the anatomical shape of the body providing excellent support for the cervical spine and ensuring its proper position during sleep.

"Notches" on both sides allow for double-sided use, both for those who require high support of 12cm, as well as for those to whom the height of 6cm is optimal. These "notches" make the pillow versatile, allowing you to choose your favourite sleeping position. Ergonomic profiling corrects the position of the spine during sleep through the support of the head, neck and shoulders.

Ultra breathable lining, made of special Mesh material ensures extraordinary comfort

Properties:

- » Due to its ergonomic profiling, pillow corrects the position of the spine during sleep through the support of the head, neck and shoulders
- » "Notches" on both sides allow for double-sided use, both for those who require high support of 12cm, as well as for those to whom the height of 6cm is optimal
- » These "notches" make the pillow versatile, allowing you to choose your favourite sleeping position
- » Relieves spinal, neck and shoulder pain by reducing pressure points
- » With thermoforming foam it perfectly adapts to the anatomical shape of the body
- » Innovative foam prevents the growth of microorganisms, allergens and absorption of dust particles
- » Ultra breathable lining, made of special Mesh material ensures extraordinary comfort
- » The pillow comes with a zipper, therefore maintaining the pillow in proper hygienic condition is not a problem

Highly recommended in the case of:

- » snoring
- » spine pain
- » shoulder pain
- » hand numbness
- » insomnia
- » nervous tensions

Pillow comes with two covers, inner and outer. Inner cover provides additional protection of the foam. Outer cover is made from innovative, ultra-breathable lining made of special Mesh material to guarantee extraordinary comfort. Zippered outer cover allows for washing of the pillow at the same time preserving its aesthetic appearance.

Dimensions: 54x33x12/6cm; **Cover:** Jersey; **Colour:** White

3 years WARRANTY

Cover available separately: cat. C-5433EP



UNIKALNY ERGONOMICZNY KSZTAŁT

Unique ergonomic design

ARmedical®
Pillow Line



MFP-6035

Profiled orthopaedic pillow "Silent Dream"

Made from an innovative thermoforming foam with a shape "memory" effect. This allows the pillow to memorise a shape of the head and as a result lower the pressure on the body keeping proper blood circulation. Thanks to the special profiles the pillow is recommended for snorers, it prevents mild symptoms of sleep apnoea, may prevent or reduce snoring.

Thanks to the ergonomic profiling, pillow corrects the position of the spine during sleep through the support of the head, neck and shoulders. Relieves spinal, neck and shoulder pain by reducing pressure points. With thermoforming foam it perfectly adapts to the anatomical shape of the body.

- » Special profiles allow to limit snoring keeping the airway open
- » Thanks to the ergonomic profiling, pillow corrects the position of the spine during sleep through the support of the head, neck and shoulders
- » Prevents morning stiffening "numbness" of hands
- » Perfectly adapts to the anatomical shape of the body

Highly recommended in the case of:

- » snoring
- » spine pain
- » shoulder pain
- » hand numbness
- » insomnia
- » nervous tensions

Pillow comes with two covers, inner and outer. Inner cover provides additional protection of the foam. Outer cover is zippered and made from innovative, super-soft micro fibre. This allows washing of the pillow at the same time preserving its aesthetic appearance.

Dimensions: 60x35x12cm; **Cover:** micro-fibre; **Colour:** Ecrú

3 years WARRANTY

Cover available separately: cat. C-6035





TFP-5030TP Traditional orthopaedic pillow "EXCLUSIVE DREAM"

Made from an innovative thermoforming foam with a shape "memory" effect. The pillow "memorises" the imprint of the head, ensuring a much lower pressure on the body without interfering blood circulation. Recommended for those who prefer traditional shape and sleep in any position.

It is a combination of regularly shaped pillow with innovative foam which alleviates the pressure.

Pillow has a super soft innovative micro-fiber pillowslip (zippered) which allows washing of the pillow at the same time preserving its aesthetic appearance.

Highly recommended in the case of:

- » spine pain
- » shoulder pain
- » hand numbness
- » snoring
- » insomnia
- » nervous tensions



Properties:

- » Perfectly adapts to the anatomical shape of the body
- » Recommended for those who prefer traditional shape and sleep in any position
- » It is a combination of regularly shaped pillow with innovative foam which alleviates the pressure
- » Due to the innovative foam with "memory" effect, the pillow corrects the position of the spine during sleep through the support of the head, neck and shoulders
- » Relieves spinal, neck and shoulder pain by reducing pressure points
- » Innovative foam prevents the growth of microorganisms, allergens and absorption of dust particles
- » Super-soft, innovative pillowslip made from micro-fibre comes with a zipper, thus maintaining the pillow in proper hygienic condition is not a problem



Dimensions: 50x30x10cm; **Cover:** Micro-fibre; **Colour:** Ecrú

3 years WARRANTY

Cover available separately: cat. C-5030TP



TFP-6040 Traditional orthopaedic pillow "EXCLUSIVE DREAM"

Made from an innovative thermoforming foam with a shape "memory" effect. This allows the pillow to memorise a shape of the head and as a result lower the pressure on the body keeping proper blood circulation.

Pillow has a super soft, innovative micro-fibre pillowslip (zippered) which allows washing of the pillow at the same time preserving its aesthetic appearance.

- » Provides optimal alignment of the cervical vertebrae during sleep
- » Thanks to the shape "memory" foam, pillow corrects the position of the spine during sleep through the support of the head, neck and shoulders
- » Perfectly adapts to the anatomical shape of the body
- » Relieves spinal, neck and shoulder pain by reducing pressure points
- » Innovative foam prevents the growth of microorganisms, allergens and absorption of dust particles

Highly recommended in the case of:

- » spine pain
- » shoulder pain
- » hand numbness
- » snoring
- » insomnia
- » nervous tensions



Dimensions: 60 x 40 x 12cm; **Colour:** Ecrú

3 years WARRANTY

Cover available separately: cat. C-6040



MFP-4026 Profiled orthopaedic pillow
"BABY DREAM"

Thanks to the unique design technology, pillow adapts perfectly to the shape of the child's body whether the child is sleeping on his/her side or on the back. The head almost submerges in the pillow's foam, thereby increasing the contact surface of the body/head with the surface as a result reducing pressure force.

Additional profiling for the head ensures proper alignment of the child in the right part of the pillow. Air holes provide constant air circulation even when the child sleeps on the side. Pillow has a super soft, innovative micro-fibre pillowslip (zippered) which allows washing of the pillow at the same time preserving its aesthetic appearance.

It provides:

- » proper sleep position from an early age
- » spine rehabilitation and its proper development
- » relief when the child is sensitive to the pressure of the neck
- » calm and healthy sleep
- » excellent cuddling


Properties:

- » Guarantees correct positioning of the cervical spine during sleep through the proper support of the head, neck and shoulders
- » Perfect for children who are more sensitive to the pressure of the neck
- » Relieves back, neck and arm pain by reducing pressure points
- » Contributes to spine rehabilitation and its proper development
- » With thermoforming foam it perfectly adapts to the anatomical shape of the body
- » Innovative foam prevents the growth of microorganisms, allergens and absorption of dust particles

Dimensions: 40x26x7cm; **Cover:** micro-fibre; **Colour:** Ecrú;

3 years WARRANTY

Cover available separately: cat. C-4026





MTP-2828 Orthopaedic pillow "EXCLUSIVE TRAVEL"

Made from an innovative thermoforming foam with a shape "memory" effect.

- » Thanks to the ergonomic profiling guarantees proper head and neck support while travelling
- » Relieves spinal, neck and shoulder pain by reducing pressure points
- » Perfectly adapts to the anatomical shape of the body
- » Innovative foam prevents the growth of microorganisms allergens and absorption of dust particles

Used for:

- » overload conditions of the cervical spine
- » instability of the cervical spine
- » degenerative spine disease
- » long trip by car, train or plane

Pillow has a super soft velur pillowslip (zippered) which allows washing of the pillow at the same time preserving its aesthetic appearance.

Dimensions: 28 x 28 x 7cm; **Colour:** Navy blue; 3 years WARRANTY



ARmedical®
Pillow Line

MTP-3230 Orthopaedic pillow "EXCLUSIVE TRAVEL"

Made from an innovative thermoforming foam with a shape "memory" effect. Thanks to that the pillow perfectly adapts to the anatomical shape of the body providing excellent support for the neck.

Thanks to the ergonomic profiling guarantees proper head and neck support while travelling. Relieves headache, spinal, neck and shoulder pain by reducing pressure points. Innovative foam prevents the growth of microorganisms, allergens and absorption of dust particles. Super-soft velour pillowcase comes with a zipper which allows to keep the pillow in good hygienic condition.

Used for:

- » overload conditions of the cervical spine
- » instability of the cervical spine
- » degenerative spine disease
- » long trip by car, train or plane

Pillow comes with two covers, inner and outer. Inner cover provides additional protection of the foam. Outer cover is zippered and made from super-soft velour. This allows washing of the pillow at the same time preserving its aesthetic appearance.

Dimensions: 32x30x12cm; **Cover:** velur, **Colour:** Navy blue



ARmedical®
Pillow Line



MFP-3433 Orthopaedic lumbar pillow "EXCLUSIVE SUPPORT"

Made from an innovative thermoforming foam with a shape "memory" effect.

- » Thanks to the ergonomic profiling, pillow corrects the position of the spine during sitting through the support of the lumbar spine
- » Helps with spine rehabilitation
- » Helps with degenerative disc disease treatment.
- » Helps relieve back pain caused by prolonged seat eg. In the office or in a car
- » Provides a pleasant wrapping feeling and "relief"
- » Additionally, it can be used as a support for the knees and legs giving the effect of relaxation for "tired" legs

With thermoforming foam it perfectly adapts to the anatomical shape of the body.

Innovative foam prevents the growth of microorganisms, allergens and absorption of dust particles.

Super-soft innovative velour pillowcase comes with a zipper which allows to keep the pillow in good hygienic condition.

Dimensions: 34x33x12cm; **Cover:** velur, navy blue colour;
 3 years WARRANTY



ARmedical®
Pillow Line

MFP-5344 Orthopaedic lumbar pillow "EXCLUSIVE SUPPORT"

Made from an innovative thermoforming foam with a shape "memory" effect. Thanks to that the pillow perfectly adapts to the anatomical shape of the body providing excellent lumbar support, ensuring its proper position during sitting.

Thanks to the ergonomic profiling, pillow corrects the position of the spine during sitting through the support of the lumbar spine. Helps with spine rehabilitation and with degenerative disc disease treatment. Helps relieve back pain caused by prolonged seat eg. In the office or in a car provides a pleasant wrapping feeling and "relief" innovative foam prevents the growth of microorganisms, allergens and absorption of dust particles. Super-soft innovative velour pillowcase comes with a zipper which allows to keep the pillow in good hygienic condition. Pillowcase comes with an elastic band for easy mounting.




Pillow provides:


- » correct position of the lumbar spine
- » reduction of pain in the lumbar area
- » relief for the spine during prolonged sitting
- » help with degenerative disc disease treatment
- » discomfort reduction associated with the overload of the spine

Pillow comes with two covers, inner and outer. Inner cover provides additional protection of the foam. Outer cover is made from innovative, super-breathable fabric which guarantees a full comfort of use.

Dimensions: 53x44x9cm; **Cover:** jersey; **Colour:** black
 3 years WARRANTY



ARmedical®
Pillow Line



MFP-2520 Orthopaedic knee "separator" pillow "EXCLUSIVE SUPPORT"

Made from an innovative thermoforming foam with a shape "memory" effect.

- » Helps to provide correct anatomical position of the hips and pelvis while lying on the side eg. In the process of rehabilitation after surgery
- » Reduces pain caused by improper arrangement of the hips during lying on the side
- » Reduces tension in the knee joints
- » Provides excellent support of the knees, thus reduces the pain in the knees caused by uncomfortable pressure. It prevents legs from sweating

With thermoforming foam it perfectly adapts to the anatomical shape of the body. Innovative foam prevents the growth of microorganisms, allergens and absorption of dust particles. Super-soft innovative velour pillowcase comes with a zipper which allows to keep the pillow in good hygienic condition.

Dimensions: 25x20x15cm; **Cover:** velur, navy blue; 3 years WARRANTY



MFP-4220 Orthopaedic lumbar pillow "EXCLUSIVE SUPPORT"

Made from an innovative thermoforming foam with a shape "memory" effect.


Highly recommended in the case of:

- » back, shoulder, neck, knee, and hips pain
- » hand numbness
- » snoring
- » insomnia
- » nervous tension, "heavy" legs syndrome

With thermoforming foam it perfectly adapts to the anatomical shape of the body. Innovative foam prevents the growth of microorganisms, allergens and absorption of dust particles. Super-soft innovative velour pillowcase comes with a zipper which allows to keep the pillow in good hygienic condition.

- » Thanks to its shape the pillow is versatile, it can be used to support neck, lumbar spine knees and between the legs as a separator
- » Helps relieve pain in the neck, back and legs – depending on how it is used

Dimensions: 42x20x10cm; **Cover:** velur, navy blue; 3 years WARRANTY



Rehabilitation balls



DSC-33 Dynamic seat cushion "DYNAPAD"

Dynamic seat cushion "DYNAPAD" is designed for correction of posture defects and for back muscles strengthening. It relieves the lumbar spine during sitting, forces correct sitting position and stimulates micro-movements of the intervertebral discs. Highly recommended for people who sit for long hours. Pressure control valve allows to adjust the hardness of the cushion. Cushion is dual-function, thanks to the bumps located on one side it can be used for feet massage in order to improve blood circulation and stimulate their receptors.

- » Perfect for correcting posture defects
- » Reduces tensions during sitting and strengthens different muscle groups
- » Perfect for feet massage and balance exercises

Diameter: 33cm

Maximum load: 200 kg

Material: PVC ECOfriendly 

PVC does not contain phthalates – laboratory tested material does not contain toxic substances and heavy metals

-  Red
-  Blue
-  Gray



FREE pump!





DSC-35 Dynamic seat cushion "DYNAPAD MAX"

Dynamic seating cushion "DYNAPAD MAX" is designed for correction of posture defects and for back muscles strengthening. It relieves the lumbar spine during sitting, forces correct sitting position and stimulates micro-movements of the intervertebral discs. Highly recommended for people who sit for long hours. Perfect for rehabilitation – during balancing exercises, training and movement coordination, and for exercises aimed at strengthening the muscles of the abdomen and legs. Pressure control valve allows to adjust the hardness of the cushion. Cushion is dual-function, thanks to the bumps located on one side it can be used for feet massage in order to improve blood circulation and stimulate their receptors. Cushion is made of anti-slip material and while doing balance exercises (standing on it) restricts the movement of the foot on its surface. The material also limits sliding of the cushions on hard ground.

- » Perfect for correcting posture defects
- » Reduces tensions during sitting and strengthens different muscle groups
- » Perfect for balance exercises
- » Perfect for feet massage

Diameter: 35cm

Height: 8cm

Maximum load: 200 kg

Material: PVC **ECOfriendly** 

PVC does not contain phthalates – laboratory tested material does not contain toxic substances and heavy metals.

FREE pump!



DSC-37 Dynamic seat cushion "DYNAPAD BALANCE"

Dynamic seat cushion "DYNAPAD BALANCE" is designed for rehabilitation, especially for the balance exercises. movement coordination training and exercises aimed at strengthening the muscles of the abdomen and legs. It is also designed for correction of posture defects and for back muscles strengthening. It relieves the lumbar spine during sitting, forces correct sitting position and stimulates micro-movements of the intervertebral discs. Highly recommended for people who sit for long hours. Pressure control valve allows to adjust the hardness of the cushion. Cushion is dual-function, thanks to the bumps located on one side it can be used for feet massage in order to improve blood circulation and stimulate their receptors.

- » Perfect for balance exercises
- » Perfect for correcting posture defects
- » Reduces tensions during sitting and strengthens different muscle groups

Diameter: 37,5cm

Maximum load: 220 kg

Material: PVC **ECOfriendly** 

PVC does not contain phthalates – laboratory tested material does not contain toxic substances and heavy metals.



FREE pump!





FREE pump!

**GM-55 to 85** Rehabilitation ball "CLASSIC"

Rehabilitation ball "Classic" is designed to be used for correction of posture defects, in therapy to improve physical condition and motor skills, to strengthen back muscles, to enhance joints mobility and to help in neurological diseases. Perfect for those who want to maintain good physical condition with all-purpose work-out.

- » Perfect for rehabilitation and balance exercises
- » Perfect for spine diseases, it corrects posture defects
- » Supports the mobility of the spine, strengthens the back muscles

Maximum load: 250 kg

Material: PVC ECOfriendly

PVC does not contain phthalates – laboratory tested material does not contain toxic substances and heavy metals.

Catalogue number	Ball size	Height of an exercising	Colour
GM-55	55cm	150 – 165 cm	
GM-65	65cm	165 – 175 cm	
GM-75	75cm	175 – 185 cm	
GM-85	85cm	over 185 cm	

ARmedical®
Active Line

ABS-55 to 85 ABS Rehabilitation ball

ABS rehabilitation ball is designed to be used for correction of posture defects, in therapy to improve physical condition and motor skills, to strengthen back muscles, to enhance joints mobility and to help in neurological diseases. Perfect for those who want to maintain good physical condition with all-purpose work-out. Recommended for pregnant women and for elderly.


- » Perfect for rehabilitation and balance exercises
- » Perfect for spine diseases, it corrects posture defects
- » Perfect for pregnant women, reduces tension in the the lower spine
- » Supports the mobility of the spine, strengthens the back muscles

Maximum load: 300 kg

Material: Anti-burst PVC ECOfriendly

ABS Rehabilitation ball is made of special material which ensures an even and slow loss of air under heavy mechanical damage to the ball, it helps to prevent tearing and sudden loss of air. Additionally, the material provides a good grip of the ball to the surface (no slips).

PVC does not contain phthalates – laboratory tested material does not contain toxic substances and heavy metals.



Catalogue number	Ball size	Height of an exercising	Colour
ABS-55	55cm	150 – 165 cm	
ABS-65	65cm	165 – 175 cm	
ABS-75	75cm	175 – 185 cm	
ABS-85	85cm	over 185 cm	



FREE pump!



ARmedical®
Active Line



BMB-55 to 75 Massage rehabilitation ball

Massage rehabilitation ball is designed to be used for correction of posture defects, in therapy to improve physical condition and motor skills, to strengthen back muscles, to enhance joints mobility and to help in neurological diseases. Perfect for those who want to maintain good physical condition with all-purpose work-out. With soft massaging bumps on the entire surface area it stimulates reflexology (receptors) zones on the skin acting as a stimulant.

- » Perfect for rehabilitation and balance exercises
- » Perfect for spine diseases, it corrects posture defects
- » Supports the mobility of the spine, strengthens the back muscles

Maximum load: 250 kg

Material: PVC ECOfriendly

PVC does not contain phthalates – laboratory tested material does not contain toxic substances and heavy metals.

Catalogue number	Ball size	Height of an exercising	Colour
BMB-55	55cm	150 – 165 cm	Orange
BMB-65	65cm	165 – 175 cm	Green
BMB-75	75cm	over 175 cm	Grey

ARmedical
Active Line



MB-5 to 10 Squeeze rehabilitation ball with spikes

Rehabilitation ball with spikes is designed for exercise, rehabilitation and hand massage. With soft bumps (spikes) it stimulates the sensory receptors, improves blood circulation. It is also used to massage different parts of the body. Thanks to its plasticity it is also used as an anti-stress ball. Comes with a valve to regulate the hardness.

- » Perfect for exercises and hand massage

Material: PVC ECOfriendly

PVC does not contain phthalates – laboratory tested material does not contain toxic substances and heavy metals.

Catalogue number	Ball size	Colour	Catalogue number	Ball size	Colour
MB-5	5cm	Orange	MB-8	8cm	Yellow
MB-6	6cm	Blue	MB-9	9cm	Red
MB-7	7cm	Green	MB-10	10cm	Grey



ARmedical
Active Line

TMB-20 Massage rehabilitation ball

Rehabilitation massage ball "MIDI" is designed for therapy to improve physical condition and movement coordination. Perfect for those who want to maintain good physical condition with all-purpose work-out. With soft massaging bumps on the entire surface area it stimulates reflexology (receptors) zones on the skin acting as a stimulant. Perfect for children to play.

- » Perfect for rehabilitation and general development exercises
- » Perfect for knees and elbows rehabilitation
- » Perfect to play

Diameter: 20cm

Maximum load: 30kg

Material: PVC transparent ECOfriendly

PVC does not contain phthalates – laboratory tested material does not contain toxic substances and heavy metals.



ARmedical
Active Line

**RLB 20-25** Rehabilitation ball MIDI REH

Rehabilitation ball "MIDI REH" is designed for therapy to improve physical condition and movement coordination, for isometric exercises which are used to counteract muscle atrophy, to maintain muscle activity of immobilised parts of the body (for injuries and fractures) and to recover the muscle growth and strength. Perfect for those who want to maintain good physical condition with all-purpose work-out. Perfect for children to play.

- » Perfect for rehabilitation and balance exercises
- » Perfect for spine diseases, it corrects posture defects
- » Supports the mobility of the spine, strengthens the back muscles

Maximum load: 30 kg

Material: PVC ECOfriendly

PVC does not contain phthalates – laboratory tested material does not contain toxic substances and heavy metals

Catalogue number	Ball size	Colour
RLB-20	20cm	
RLB-25	25cm	

FREE pump!



ARmedical®
Active Line

PSB-50 ABS Rehabilitation ball "PEANUT"

Rehabilitation ball "PEANUT" is designed to be used for correction of posture defects, in therapy to improve physical condition and motor skills, to strengthen back muscles, to enhance joints mobility and to help in neurological diseases. Perfect for those who want to maintain good physical condition with all-purpose work-out.

- » Perfect for rehabilitation and balance exercises
- » Perfect for spine diseases, it corrects posture defects
- » Supports the mobility of the spine, strengthens the back muscles

Maximum load: 350 kg

Size: 50x100cm

Material: Anti-burst PVC ECOfriendly

Rehabilitation ball "PEANUT" is made of special material which ensures an even and slow loss of air under heavy mechanical damage to the ball, it helps to prevent tearing and sudden loss of air. Additionally, the material provides a good grip of the ball to the surface (no slips).

PVC does not contain phthalates – laboratory tested material does not contain toxic substances and heavy metals.



FREE pump!

ARmedical®
Active Line



HB2-45 to 55 Rehabilitation ball "HOPPER" with horn handles


Rehabilitation ball "HOPPER" for children is designed to be used for correction of posture defects, in therapy to improve physical condition and motor skills, to strengthen back muscles, to enhance joints mobility and to help in neurological diseases. Perfect to improve physical condition with all-purpose work-out. Perfect to play with in order to correct posture defects in children.

- » Perfect for rehabilitation and balance exercises
- » Perfect for spine diseases, it corrects posture defects
- » Supports the mobility of the spine, strengthens the back muscles
- » Perfect to play

Maximum load: 80 kg

Material: PVC ECOfriendly

PVC does not contain phthalates – laboratory tested material does not contain toxic substances and heavy metals.



FREE pump!



Catalogue number	Ball size	Height of an exercising	Colour
HB2-45	45cm	130 – 140 cm	Red
HB2-50	50cm	140 – 150 cm	Yellow
HB2-55	55cm	150 – 165 cm	Green

ARmedical
Active Line

HB1-45 to 55 Rehabilitation ball "HOPPER" with handle

Rehabilitation ball "HOPPER" for children is designed to be used for correction of posture defects, in therapy to improve physical condition and motor skills, to strengthen back muscles, to enhance joints mobility and to help in neurological diseases. Perfect to improve physical condition with all-purpose work-out. Perfect to play with in order to correct posture defects in children.

- » Perfect for rehabilitation and balance exercises
- » Perfect for spine diseases, it corrects posture defects
- » Supports the mobility of the spine, strengthens the back muscles
- » Perfect to play

Maximum load: 80 kg

Material: PVC ECOfriendly

PVC does not contain phthalates – laboratory tested material does not contain toxic substances and heavy metals.



FREE pump!



Catalogue number	Ball size	Height of an exercising	Colour
HB1-45	45cm	130 – 140 cm	Red
HB1-50	50cm	140 – 150 cm	Yellow
HB1-55	55cm	150 – 165 cm	Green

ARmedical
Active Line

Rehabilitation and assistive equipment

AR-920 Anti Decubitis Air Bubble Mattress & Pump System

It is used as an overlay to a typical mattress on the patient's bed, and afterwards it is covered with a sheet. Two separate sections of the mattress are alternately filled by a special pump.

Thanks to this alternating operation (mattress + pump) the individual areas of the patient's body are not in a continuous contact with the ground for longer than a few minutes. The air pressure in the mattress is adjusted (increased or decreased) depending on the weight of the patient.

- » **Low noise** is a great advantage of ARmedical pump. Regardless of the load, the pump runs silently.
- » An additional advantage of this system (pump + mattress) is failure-free functioning. The package includes a repair kit (patches + glue) used in the case of a cut with sharp tool.
- » Recommended for bedridden patients to prevent bed sores and to treat existing ones.

Specifications:

- » pump with adjustable pressure
- » dimensions when inflated: 200cm x 90cm x 7 cm (plus two special extensions to put under the mattress)
- » compressor working pressure: 40 do 100mmHg (smooth adjustment)
- » full cycle time: 6 minutes
- » special hooks used for hanging the pump on the bed

Mattress available separately: cat. AR-921

Maximum load: 135kg

New quality



New quality

AR-950 Anti Decubitis Air Tubing Mattress & Pump System

It is used as an overlay to a typical mattress on the patient's bed, and afterwards it is covered with a sheet. Chambers (tubes) of the mattress are alternately filled by a special pump.

Thanks to this alternating operation (mattress + pump) the individual areas of the patient's body are not in a continuous contact with the ground for longer than a few minutes. The air pressure in the mattress is adjusted (increased or decreased) depending on the weight of the patient.

- » **Low noise** is a great advantage of ARmedical pump. Regardless of the load, the pump runs silently.
- » An additional advantage of this system (pump + mattress) is failure-free functioning. Recommended for bedridden patients to prevent bed sores and to treat existing ones.

Pump specifications:

- » type: pressure-swing
- » smooth pressure adjustment
- » working pressure: from 40 to 100mmHg
- » full cycle time: 6 minutes
- » special hooks used for hanging on the bed

Cover specifications:

- » material: nylon + PVC
- » chambers: 20
- » dimensions when inflated: 200cm x 86cm x 10cm
- » Mounting: on an existing mattress (2 overlaps)

Mattress available separately: cat. AR-951

Maximum load: 140kg



AR-018 Pedal exerciser

Pedal exerciser comes with anti-slip caps which make it possible to set it up anywhere with no need for additional racks. It is equipped with very simple and smooth resistance adjustment which works thanks to bolt mechanism. It is an ideal rehabilitation device for people with upper and lower limbs disabilities.

Use:

AR-018 Rehabilitation pedal exerciser is a manual device designed for rehabilitation. Recommended for strengthening of upper and lower limbs, especially in the case of paresis caused by stroke, accident and other various neurological disorders. The use of pedal exerciser helps to build up muscle strength, improve blood circulation in the limbs, stimulate the nervous system and improve joint mobility. Highly recommended for exercises in the case of paralysis of the peripheral nerves of the lower limbs. Use ensures:

- » efficiency improvement of cardio-respiratory system
- » blood circulation improvement in the limbs (arterial, venous and lymphatic circulation)
- » prevents muscular atrophy
- » prevents arterial and venous blood clots
- » prevents contractures in the joints of the limbs
- » exercises of the lower limbs can be performed while seated on a wheelchair or any other chair
- » exercises of the upper limbs can be performed after setting the pedal exerciser on the table

Technical specification:

- » Frame: steel powder coated
- » Pedals: Nylon
- » Anti-slip pads ABS
- » Pedalling resistance: smooth adjustment
- » Dimensions: 49x41x30cm
- » Weight: 2kg
- » Maximum total load: 70kg
- » Commercial package: 1pc
- » Collective carton box: 8pcs



Resistance adjustment



Feet straps



Available colours:



blue



silver



black

ARmedical
HomeCare Line

AR-019 Pedal exerciser – folding with meter

Pedal exerciser comes with anti-slip caps which make it possible to set it up anywhere with no need for additional racks. It is equipped with very simple and smooth resistance adjustment which works thanks to bolt mechanism. Multifunction display allows to monitor workout time, the amount of burnt calories, and cadence. It is an ideal rehabilitation device for people with upper and lower limbs disabilities.

Folding option allows for easy storage of the rotor after completed exercises.

Use:

AR-019 Rehabilitation pedal exerciser is a manual device designed for rehabilitation. Recommended for strengthening of upper and lower limbs, especially in the case of paresis caused by stroke, accident and other various neurological disorders. The use of rotor helps to build up muscle strength, improve blood circulation in the limbs, stimulate the nervous system and improve joint mobility. Highly recommended for exercises in the case of paralysis of the peripheral nerves of the lower limbs. Use ensures:

- » efficiency improvement of cardio-respiratory system
- » blood circulation improvement in the limbs (arterial, venous and lymphatic circulation)
- » prevents muscular atrophy
- » prevents arterial and venous blood clots
- » prevents contractures in the joints of the limbs
- » exercises of the lower limbs can be performed while seated on a wheelchair or any other chair
- » exercises of the upper limbs can be performed after setting the rotor on the table



Folding



Electronic display



Available colours:



blue



silver



black

Technical specification:

- » Frame: steel powder coated, folding
- » Pedals: Nylon
- » Anti-slip pads ABS
- » Pedalling resistance: smooth adjustment
- » Accessories: multifunction display
- » Dimensions: 50x36x26cm
- » Weight: 2.3kg
- » Maximum total load: 100kg
- » Commercial package: 1pc.
- » Collective carton box: 4pcs.

ARmedical
HomeCare Line

AR-010 Elbow crutch – aluminium

Elbow crutches are designed for people who need support while learning to walk after injuries of the lower limbs, to maintain the balance and weight transfer in orthopaedic trauma disorders (fractures, dislocations, sprains). They are made of two durable and lightweight aluminium tubes, plastic grip and a rubber tip (anti-slip). By using properly soft rubber compound, non-slip tip perfectly absorbs shock. Powder coating technology ensures varnish durability and aesthetic appearance for a long time.

Technical specs:

- » durable aluminium tubes
- » comfortable profiled grip
- » firm opened cuff
- » high quality rubber tip
- » cutting height adjustment: 2.5 cm
- » adjustment range: from 70.5 to 93.5 cm
- » Weight: 0.50 kg



AR-011 – rubber tips available separately:



AR-013 Aluminium T-shape cane


Walking stick is a device used by people who require support while learning to walk after suffering from lower limbs injuries, as well as to maintain the balance and to transfer body weight in the case of orthopaedic trauma and disorders. Used frequently by the elderly, the disabled and people with problems of musculoskeletal system.

It is made of durable and lightweight aluminium tubes, ergonomically contoured handle and a non-slip soft rubber tip. Non-slip soft rubber tip perfectly absorbs shock thanks to the appropriate soft rubber compound.

Painting/anodizing technology ensures durability of the paint and aesthetic look for a long time.

Technical specification:

- » durable aluminum tubes
- » comfortable, ergonomically shaped handle
- » universal handle (right or left hand)
- » additional „string” to keep the stick up
- » high quality soft rubber tip
- » adjustable every 2.5 cm
- » adjustment range from 71.0 to 94.0 cm
- » Weight: 0.30 kg
- » Weight capacity: 100kg



Available colours: Silver Copper



AR-014 Aluminium T-shape cane with soft handgrip


Walking stick is a device used by people who require support while learning to walk after suffering from lower limbs injuries, as well as to maintain the balance and to transfer body weight in the case of orthopaedic trauma and disorders. Used frequently by the elderly, the disabled and people with problems of musculoskeletal system.

It is made of durable and lightweight aluminium tubes, ergonomically contoured soft-handgrip and a non-slip rubber tip. Non-slip rubber tip perfectly absorbs shock thanks to the appropriate soft rubber compound.

Powder coating technology ensures durability of the paint and aesthetic look for a long time.

Technical specification:

- » durable aluminum tubes
- » comfortable, soft, ergonomically shaped handgrip
- » soft, universal handle (right or left hand)
- » additional „string” to keep the stick up
- » high quality soft rubber tip
- » adjustable height every 2.5 cm
- » range adjustment from 69.0 to 92.0 cm
- » Weight: 0.37 kg
- » Weight capacity: 100kg



Available colours: Anthracite



AR-015 Aluminium T-shape cane – foldable

Walking stick is a device used by people who require support while learning to walk after suffering from lower limbs injuries, as well as to maintain the balance and to transfer body weight in the case of orthopaedic trauma and disorders. Used frequently by the elderly, the disabled and people with problems of musculoskeletal system.

It is made of durable and lightweight aluminium tubes, ergonomically contoured handle and a non-slip rubber tip. Non-slip ferrule perfectly absorbs shock thanks to the appropriate soft rubber compound. Painting/anodizing technology ensures durability of the paint and aesthetic look for a long time. With folding option and low weight, the stick can be easily put in a bag or backpack. It will always be at hand when needed.

Technical specification:

- » durable aluminium tubes (4 parts assembly)
- » comfortable, ergonomically shaped handle
- » universal handle (right or left hand)
- » additional „string” to keep the stick up
- » high quality soft rubber ferrule
- » adjustable height every 2.5 cm
- » range adjustment from 82.0 to 93.0 cm
- » Weight: 0.37 kg
- » Weight capacity: 100kg




Available colours: Silver



Copper



**AR-016** Aluminium T-shape Tripod Cane

Walking stick is a device used by people who require support while learning to walk after suffering from lower limbs injuries, as well as to maintain the balance and to transfer body weight in the case of orthopaedic trauma and disorders. Used frequently by the elderly, the disabled and people with problems of musculoskeletal system

It is made of durable and lightweight aluminium tubes with a three-point support, ergonomically contoured handle and a non-slip rubber tips. Non-slip rubber tips perfectly absorbs shock thanks to the appropriate soft rubber compound. Painting/anodizing technology ensures durability of the paint and aesthetic look for a long time.

Thanks to the wide base it is designed for people requiring three-point, solid support.

Technical specification:

- » durable aluminium tubes, chromed steel base
- » comfortable, ergonomically shaped handle
- » universal handle (right or left hand)
- » high quality soft rubber tips
- » adjustable height every 2.5 cm
- » range adjustment from 74.0 to 97.0 cm
- » Weight: 0.7 kg
- » Weight capacity: 100kg



Available colour: Silver



ARmedical®
HomeCare Line

AR-017 Aluminium T-shape Quad Cushion Cane with soft handgrip

Walking stick is a device used by people who require support while learning to walk after suffering from lower limbs injuries, as well as to maintain the balance and to transfer body weight in the case of orthopaedic trauma and disorders. Used frequently by the elderly, the disabled and people with problems of musculoskeletal system.

It is made of durable and lightweight aluminium tubes. The base is made of extremely durable material with a four-point support, soft and ergonomically contoured handgrip. The base is equipped with anti-slip rubber tips. Non-slip rubber tips perfectly absorb shock thanks to the appropriate soft rubber compound. Painting/anodizing technology ensures durability of the paint and aesthetic look for a long time. Thanks to the wide base it is designed for people requiring four-point, solid support.

Technical specification:

- » durable aluminum tubes
- » comfortable, soft, ergonomically shaped handgrip
- » soft, universal handle (right or left hand)
- » additional „string“ to keep the stick up
- » durable base with soft rubber tips
- » adjustable height every 2.5 cm
- » range adjustment from 62.0 to 85.0 cm
- » Weight: 0.45 kg
- » Weight capacity: 100kg



ARmedical®
HomeCare Line

Available colours: Silver



Anthracite






AR-017W Aluminium T-shape High Quad Cushion Cane

Walking stick is a device used by people who require support while learning to walk after suffering from lower limbs injuries, as well as to maintain the balance and to transfer body weight in the case of orthopaedic trauma and disorders. Used frequently by the elderly, the disabled and people with problems of musculoskeletal system.

It is made of durable and lightweight aluminium tubes with a four-point support, ergonomically contoured handle and a non-slip rubber tips. Non-slip rubber tips perfectly absorb shock thanks to the appropriate soft rubber compound. Painting/anodizing technology ensures durability of the paint and aesthetic look for a long time. Thanks to the wide base it is designed for people requiring four-point, very solid support.

Technical specification:

- » durable aluminium tubes, chromed steel base
- » comfortable, ergonomically shaped handgrip
- » universal handle (right or left hand)
- » high quality soft rubber tips
- » adjustable height every 2.5 cm
- » range adjustment from 74.0 to 97.0 cm
- » Weight: 0.8 kg
- » Weight capacity: 100kg



Available colours: Silver 



AR-001 Aluminium folding walker. Rolling – stable. Single H-shaped frame

General information:

Rollator is a physical-rehab medical device for disabled people with musculoskeletal system disorder, with verticalisation problems or having difficulties maintaining balance while walking. It is also used to maintain balance while learning to walk after a serious illness or injury.

Product information:

- » reinforced aluminium frame
- » folded with one button
- » dual-function: rolling-stable
- » stable and robust design
- » adjustable height
- » ergonomic grips

Technical specs:

- » Width: 50 cm
- » Min. height: 79 cm
- » Max. height: 96 cm
- » Weight: 2.3 kg
- » Maximum load: 136 kg

Available colours: Silver  Copper 



AR-002**Aluminium folding walker Rolling – stable.
Double H-shaped frame****General information:**

Rollator is a physical-rehab medical device for disabled people with musculoskeletal system disorder, with verticalisation problems or having difficulties maintaining balance while walking. It is also used to maintain balance while learning to walk after a serious illness or injury.

Product information:

- » double-reinforced aluminium frame
- » folded with one button
- » dual-function: rolling-stable
- » stable and robust design
- » adjustable height
- » ergonomic grips

Technical specs:

- » Width: 50 cm
- » Min. height: 79 cm
- » Max. height: 96 cm
- » Weight: 2.5 kg
- » Maximum load: 136 kg



ARmedical®
HomeCare Line

Available colours:**AR-003****Aluminium folding walker
with two 5-inch wheels****General information:**

Rollator is a physical-rehab medical device for disabled people with musculoskeletal system disorder, with verticalisation problems or having difficulties maintaining balance while walking. It is also used to maintain balance while learning to walk after a serious illness or injury.

Product information:

- » reinforced aluminium frame
- » folded with one button
- » stable and robust design
- » adjustable height
- » ergonomic grips
- » front wheels 5" (12 cm)

Technical specs:

- » Width: 50 cm
- » Min. height: 70 cm
- » Max. height: 83 cm
- » Weight: 2.3 kg
- » Maximum load: 136 kg



ARmedical®
HomeCare Line

Available colours:



AR-008 Aluminium folding walker – 3 functional

General information:

Walker is a physical-rehab medical device for disabled people with musculoskeletal system disorder, with verticalisation problems or having difficulties maintaining balance while walking. It is also used to maintain balance while learning to walk after a serious illness or injury.

Product information:

- » reinforced aluminium frame
- » folded with one button
- » stable and robust design
- » adjustable height
- » ergonomic hand grips
- » front wheels 5" (12 cm) + additional footer to stepping-stable function

Technical specs:

- » Width: 50 cm
- » Min. height: 70 cm
- » Max. height: 83 cm
- » Weight: 2.3 kg
- » Maximum load: 136 kg



Available colours:



AR-009 Aluminium folding walker with ball – Two-wheels

General information:

Walker is a physical-rehab medical device for disabled people with musculoskeletal system disorder, with verticalisation problems or having difficulties maintaining balance while walking. It is also used to maintain balance while learning to walk after a serious illness or injury.

Product information:

- » aluminium frame
- » fold by pulling the ball at the front of the frame
- » stable and robust design
- » adjustable height
- » ergonomic handles
- » front wheels 5" (12 cm)

Technical specification:

- » Width: 60 cm
- » Min. height: 78 cm
- » Max. height: 85 cm
- » Weight: 1.9 kg
- » Weight capacity: 110 kg



Available colour:



AR-020 Two-wheel rehabilitation walker STANDARD

Walker is a physical-rehab medical device for disabled people with musculoskeletal system disorder, with verticalisation problems or having difficulties maintaining balance while walking. It is also used to maintain balance while learning to walk after a serious illness or injury. Used frequently by the elderly, the disabled and people with problems of musculoskeletal system.

The great advantages of this walker are its low weight and small dimensions when folded, making it easy to store. Equipped with a seat to increase the comfort of the user - a chance to rest during long walks.

Product information:

- » powder coated, folding aluminium frame
- » stable and durable design
- » with height-adjustable handles (every 2.5cm)
- » contoured handles
- » 2 maintenance-free wheels, full – 10 cm
- » 2 legs with non-slip rubber tips
- » soft seat, covered with PVC

Technical specification:

- » Seat dimensions: 42 x 21cm
- » Height to the seat: 52cm
- » Internal width between the tubes: 43.5cm
- » Handles height: 75,5-92cm (every 2,5cm)
- » Overall width: 49cm
- » Overall length: 63cm
- » Width (thickness when folded): 18cm
- » Weight: 4.5 kg
- » Weight capacity: 136kg

Available colour: Blue



ARmedical®
HomeCare Line

**AR-021** Two-wheel walker DELUXE

Walker is a physical-rehab medical device for disabled people with musculoskeletal system disorder, with verticalisation problems or having difficulties maintaining balance while walking. It is also used to maintain balance while learning to walk after a serious illness or injury. Used frequently by the elderly, the disabled and people with problems of musculoskeletal system.

The great advantages of this walker are its low weight and small dimensions when folded, making it easy to store. Equipped with a seat and a backrest to increase the comfort of the user - a chance to rest during long walks. Dual front wheels ensure greater stability.

Product information:

- » foldable alu-steel frame
- » powder coated
- » modern design
- » stable and durable design
- » with height-adjustable handles (every 2.5cm)
- » ergonomically shaped handles
- » 2 front, maintenance-free wheels, full – 10 cm
- » 2 legs with non-slip rubber tips
- » equipped with seat and backrest
- » seat base is made of plastic

Technical specification:

- » Seat dimensions: 37 x 24cm
- » Height to the seat: 49cm
- » Internal width between the tubes: 45cm
- » Handles height: 79-88cm (every 2,5cm)
- » Overall width: 59cm
- » Overall length: 61cm
- » Width (thickness when folded): 21cm
- » Weight: 4.8 kg
- » Weight capacity: 136kg

Available colour: Graphite



ARmedical®
HomeCare Line



AR-004 3-wheel Rollator with storage bag and basket

General information:

Rollator is a physical-rehab medical device for disabled people with musculoskeletal system disorder, with verticalisation problems or having difficulties maintaining balance while walking. It is also used to maintain balance while learning to walk after a serious illness or injury.

Product information:

- » powder coated, foldable aluminium frame
- » stable and robust design
- » with height-adjustable handles
- » ergonomically sculptured handles
- » front maintenance-free wheels 8" (20 cm)
- » standard equipment:
 - parking brake
 - storage bag with zipper
 - detachable basket with cover

Technical specs:

- » Min. height: 82 cm
- » Max. height: 98 cm
- » Weight: 8.3 kg
- » Maximum load: 110 kg



Available colour: Navy blue



AR-005 4-wheel steel rehab Rollator CLASSIC

General information:

Rollator is a physical-rehab medical device for disabled people with musculoskeletal system disorder, with verticalisation problems or having difficulties maintaining balance while walking. It is also used to maintain balance while learning to walk after a serious illness or injury.

Product information:

- » powder coated, foldable aluminium frame
- » stable and robust design
- » with height-adjustable handles
- » ergonomically sculptured handles
- » front maintenance-free wheels 8" (20 cm)
- » standard equipment:
 - parking brake
 - front shopping basket
 - tray mounted on the seat
 - plastic seat

Technical specs:

- » Seat width: 38 cm
- » Seat height: 61 cm
- » Grips height: 75 cm
- » Max. grips height: 98 cm
- » Weight: 11.9 kg
- » Maximum load: 136 kg



Available colour: Blue



FEATURE
PARKING BRAKE

foldable

AR-006**4-wheel aluminium rehab Rollator
PRESTIGE****General information:**

Rollator is a physical-rehab medical device for disabled people with musculoskeletal system disorder, with verticalisation problems or having difficulties maintaining balance while walking. It is also used to maintain balance while learning to walk after a serious illness or injury.

Product information:

- » powder coated, foldable aluminium frame
- » stable and robust design
- » with height-adjustable handles
- » ergonomically sculptured handles
- » 4 maintenance-free foam wheels, full – 8" (20 cm)
- » standard equipment:
 - parking brake
 - front shopping basket
 - tilted plastic seat
 - foldable and detachable backrest

Technical specs:

- » Seat width: 37 cm
- » Min. seat height: 52 cm
- » Max. seat height: 60 cm
- » Min. handles height: 94 cm
- » Max. handles height: 101 cm
- » Weight: 6 kg
- » Maximum load: 136 kg

max
136
kg**Available colour:** Maroon**AR-007****4-wheel aluminium rehab Rollator DYNAMIC****General information:**

Rollator is a physical-rehab medical device for disabled people with musculoskeletal system disorder, with verticalisation problems or having difficulties maintaining balance while walking. It is also used to maintain balance while learning to walk after a serious illness or injury.

Easy transport and storage**Product information:**

- » powder coated, foldable aluminium frame
- » ultrafast folding with the strap located on the seat
- » stable and robust design
- » with height-adjustable handles
- » ergonomically sculptured handles
- » 4 maintenance-free foam wheels, full – 8" (20 cm)
- » standard equipment:
 - parking brake
 - detachable vinyl storage bag with shoulder pad
 - vinyl seat
 - vinyl backrest, detachable

Parametry techniczne:

- » Seat width: 45 cm
- » Seat height: 56 cm
- » Min. handles height: 81 cm
- » Max. handles height: 94 cm
- » Weight: 7 kg
- » Maximum load: 136 kg

max
136
kg**Available colour:** Silver nickel

foldable

PARKING
BRAKE FEATURE



AR-022

Four-wheel aluminium rollator with footrest COMFORT

Walker is a physical-rehab medical device for disabled people with musculoskeletal system disorder, with verticalisation problems or having difficulties maintaining balance while walking. It is also used to maintain balance while learning to walk after a serious illness or injury. Used frequently by the elderly, the disabled and people with problems of musculoskeletal system.

The great advantage of this walker is the transport function. Detachable backrest, seat and folding footrest allow for patient transportation. Moreover, the seat with backrest significantly improve the comfort - a chance to rest during long walks. In addition, folding option makes it easy to store.

Product information:

- » coated aluminium frame - foldable
- » stable and durable design
- » with height-adjustable handles (every 2.5cm)
- » contoured handles
- » 4 maintenance-free wheels, full – 20 cm
- » swivel front wheels
- » standard equipment:
 - parking brake
 - shopping basket
 - plastic seat
 - soft, detachable backrest
 - folding footrest

Technical specification:

- » Seat dimensions: 35.5 x 33cm
- » Height to the seat: 54cm
- » internal width between the tubes: 46cm
- » Handles height: 79-91cm
- » Overall width: 60cm
- » Overall length (folded footrest): 67.5cm
- » Overall length (unfolded footrest): 92cm
- » Width (thickness when folded): 33cm
- » Weight: 8.2 kg
- » Weight capacity: 136kg

Available colour: Maroon



PARKING BRAKE



DETACHABLE BACKREST



FOLDING FOOTREST



Ergonomic



Brake for the caregiver with lock function



ultralight
13kg

AR-300

**Aluminium wheelchair
ERGONOMIC**

The wheelchair is a medical device intended for people with disabilities. Ergonomic is an ultralight aluminium wheelchair made using advanced technology. The wheelchair can be easily folded thanks to its lightweight, x-shaped design and "folding backrest". Another convenience is the "quick-release" system used in rear wheels which allows to quickly disassemble them (one-button). A special anti-bacterial, removable and washable material on the seat and backrest is extremely practical in everyday use. With the accessory and robust standard of craftsmanship the wheelchair will be very practical and useful in everyday life of a disabled person or a person having a problem with a body part movement.

Product information:

- » modern design
- » folding, ultra-lightweight aluminium x-shaped frame
- » folding backrest
- » soft, antibacterial, removable, washable seat and backrest
- » swing-away and removable footrests with adjustable height
- » footrests equipped with safety straps
- » folding and rotated armrests
- » hip belt to protect a person from falling off the wheelchair
- » maintenance-free front wheels
- » maintenance-free rear wheels - quick-release (one-button disassemble system)
- » aluminium handrims
- » dual braking system (for the patient and caregiver)
- » brake for the caregiver with lock function
- » brake for the patient with pressure regulation
- » footrest for the caregiver to assist with foot-lifting of the wheelchair (to pass over a kerb)
- » frame colour: silver metallic

Technical specification:

- » Weight capacity: 114kg
- » Seat width: 41 cm or 46cm
- » Weight: 13kg
- » Front wheels: maintenance-free: 18x 4,5cm (7" x 1-3/4")
- » Rear wheels: maintenance-free: 61 x 2,5cm (24"x1")
- » Armrests: folding and rotated
- » Footrests: detachable, swing-away
- » Backrest: Folding
- » Total external dimensions: 107x68x92cm



Soft, antibacterial, detachable, washable seat and backrest



Folding armrest




Quick-release rear wheels



Swing-away and detachable footrests



Brake for the patient with pressure regulation



	Seat size [cm]							Footrests	Armrests	
	D	A	B	C	E	F	G	H	J	K
Width 16" (41cm)	41	107	63	92	38-43	48	48	38-51	23	30
Width 18" (46cm)	46	107	68	92	38-43	48	48	38-51	23	30



AR-330 Aluminium wheelchair DYNAMIC

The wheelchair is a medical device intended for people with disabilities. Dynamic is an ultralight aluminium wheelchair made using advanced technology. The wheelchair can be easily folded and transported in a car thanks to its lightweight, x-shaped design and "quick-release" rear wheels. "Quick release" is a quick disassemble system (one-button). Soft nylon seat, backrest and additional cushion ensure both comfort as well as the ability to keep the wheelchair in a proper hygienic condition. With the accessory (safety belt, brake for the patient and caregiver, detachable seat cushion) and robust standard of craftsmanship, the wheelchair will be very practical and useful in everyday life of a disabled person or a person having a problem with a body part movement.

Product information:


- » modern design
- » folding, ultra-lightweight, aluminium, double x-shaped frame
- » soft, nylon seat and backrest (bag/pocket on the back)
- » soft seat cushion
- » swing-away and removable footrests with adjustable height
- » footrests equipped with safety straps
- » folding armrests
- » hip belt to protect a person from falling off the wheelchair
- » maintenance-free front wheels
- » pneumatic rear wheels with "quick-release" (one-button disassemble system)
- » aluminium handrims
- » dual braking system (for the patient and caregiver)
- » footrest for the caregiver to assist with foot-lifting of the wheelchair (to pass over a kerb)
- » ergonomic handles for caregiver
- » reflectors on the rear wheels
- » frame colour: Silver metallic

Technical specification:

- » Weight capacity: 130kg
- » Seat width: 41cm, 46cm, 50cm
- » Weight: 16kg
- » Front wheels: maintenance-free: 20x5cm (8"x 2")
- » Rear wheels: pneumatic: 61x2,5cm (24"x1")
- » Armrests: folding, PU material
- » Footrests: detachable, swing-away with adjustable height
- » Total external dimensions: 106x62x93cm



Brake for the caregiver



Pocket on the back



Brake for the patient



Quick-release rear wheels






Folding armrest



Swing-away and detachable footrests



Footrests with height adjustment



Optional: Anti-tipping wheels
Num cat.: AR-330-ANTY



	Seat size [cm]							Footrests	Armrests	
	D	A	B	C	E	F	G	H	J	K
Width 16" (41cm)	41	106	58	93	44	44	54	30-41	37	30
Width 18" (46cm)	46	106	62	93	44	44	54	30-41	37	31
Width 20" (50cm)	50	106	66	93	44	44	54	30-41	37	32

max **135** kg
hardy

AR-350

**Aluminium wheelchair
 PRESTIGE**



Quick-release rear wheels



Armrests height adjustment



Handles height adjustment



Footrests with height adjustment



Anti-flip wheels with height adjustment






Mini spare wheels



Rear wheels height adjustment



Front wheels height adjustment



The wheelchair is a medical device intended for people with disabilities. Prestige is a lightweight, modern aluminium wheelchair made using highly advanced technologies, possessing an extremely practical and functional standard equipment. The wheelchair can be easily folded and transported in a car thanks to its lightweight, x-shaped design and "quick-release" rear wheels. "Quick release" is a quick disassemble system (one-button). Soft nylon seat, backrest and additional cushion ensure both comfort as well as the ability to keep the wheelchair in a proper hygienic condition. With the accessory (adjustable centre of gravity, adjustable handles for caregiver, folding armrests with adjustable height, mini spare wheels, anti-tip wheels with height adjustment, safety belt, brake for the patient, detachable seat cushion) and robust standard of craftsmanship, the wheelchair will be very practical and useful in everyday life of a disabled person or a person having a problem with a body part movement.

Product information:

- » modern design, attractive look
- » folding, ultra-lightweight aluminium x-shaped frame
- » soft, nylon seat and backrest
- » soft seat cushion
- » hip belt to protect a person from falling off the wheelchair
- » adjustable centre of gravity (the possibility to use by the amputee)
- » height-adjustable front and rear wheels
- » folding armrests with height adjustment
- » height adjustment handles for the caregiver
- » ergonomic handles for caregiver
- » anti-tip wheels with adjustable height
- » mini spare wheels on the back
- » brake for the patient
- » swing-away and removable footrests with adjustable height
- » footrests equipped with safety straps
- » footrests equipped with an additional belt for calf support
- » maintenance-free front wheels
- » pneumatic rear wheels with "quick-release" (one-button disassemble system)
- » reflectors on the rear wheels
- » air-pump included
- » frame colour: Silver metallic

Technical specification:

- » Weight capacity: 135kg
- » Seat width: 41cm, 43cm, 46cm, 48cm, 50cm
- » Weight: 16.5kg
- » Front wheels: maintenance-free: 18x4.5cm (7" x 1-3/4")
- » Rear wheels: pneumatic: dimensions: 61x2,5cm (24" x 1")
- » Armrests: folding with height adjustment, PU material
- » Footrests: detachable, swing-away with adjustable height
- » Total external dimensions: 104x60-69x96-103cm

max **135** kg



	Seat size [cm]							Footrests	Armrests	
	D	A	B	C	E	F	G	H	J	K
Width 16" (41cm)	41	104	60	96-103	42	45-52	53	34-42	24-32	30
Width 17" (43cm)	43	104	62	96-103	42	45-52	53	34-42	24-32	31
Width 18" (46cm)	46	104	65	96-103	42	45-52	53	34-42	24-32	32
Width 19" (48cm)	48	104	67	96-103	42	45-52	53	34-42	24-32	33
Width 20" (50cm)	50	104	69	96-103	42	45-52	53	34-42	24-32	34



AR-101 Folding – steel toilet chair

General information:

Bedside Commode seat is designed for disabled people with musculoskeletal system disorder, for elderly and for those who have a problem with the use of a standard toilet bowl. Facilitates and assists during daily toilet routine. Bedside Commode can also be used as safety rails to assisting during standing up from a standard toilet bowl.

Product information:


- » powder coated, steel frame
- » foldable frame (easy storage)
- » stable and robust design
- » adjustable height
- » ergonomic armrests

Technical specs:

- » Seat width: 50 cm
- » Seat depth: 37cm
- » Toilet seat width: 35 cm
- » Toilet seat depth: 37 cm
- » Bucket capacity: 7 litres
- » Min. height: 40 cm
- » Max. height: 55 cm
- » Weight: 6 kg
- » Maximum load: 136 kg



Available colour: Gray



AR-110 / AR-115 Raised toilet seat with lid

General information:


Raised toilet seat with lid is designed for disabled people with musculoskeletal system disorder, especially for those who have a problem with standing up and sitting down. For elderly, for people with arthritis, after hip surgery. and for those who have a problem with the use of a regular-height toilet bowl. Facilitates and assists during daily toilet routine.

Attaches securely to most regular shaped toilets with easy tool-free installation. The seat is moulded from profiled plastic which is easy to clean and disinfect. Thanks to its dimensional standard it can be attached securely to most regular shaped toilets. Adjustable side clamps provide a stable placement of the seat on the toilet.

Front profile provides comfort for both women and men. Includes a lid for covering the toilet when it is not used. Ergonomic design provides comfort of use.

Technical specs

Spec:	Model AR-110	Model AR-115
Overall width:	36,5 cm	36,5 cm
Overall length:	40 cm	40 cm
Opening width:	27 cm	27 cm
Opening length:	30,5 cm	30,5 cm
Height:	10 cm	15 cm



Available colour: white



**AR-201** Rounded shower stool**General information:**

Shower Stool as an assistive device is designed for the disabled with musculoskeletal system disorder, for elderly, as well as those for which using a standard bath or shower is problematic. Helps and assists during daily toilet activities eg. using washbasin, tub or shower stall.

Product information:

- » anodized aluminium frame
- » seat made of an easily washable polyethylene
- » non-slip rubber tips
- » stable and robust design

Technical specs:

- » Seat diameter: 33 cm
- » Min. height: 38 cm
- » Max. height: 56 cm
- » Weight: 1.5 kg
- » Maximum load: 110 kg

Available colour: White

**AR-202** Shower chair without backrest**General information:**

Shower chair as an assistive device is designed for the disabled with musculoskeletal system disorder, for elderly, as well as those for which using a standard bath or shower is problematic. Helps and assists during daily toilet activities eg. using washbasin, tub or shower stall.

Product information:

- » anodized aluminium frame
- » seat made of an easily washable polyethylene
- » equipped with handles
- » non-slip rubber tips
- » stable and robust design
- » adjustable height

Technical specs:

- » Seat dimensions: 50 x 30 cm
- » Min. height: 36 cm
- » Max. height: 47 cm
- » Weight: 1.8 kg
- » Maximum load: 110 kg

Available colour: White

**AR-209** Shower stool with U-cut**General information:**

Shower stool as a medical aid is designed for the disabled with locomotor disorders, the elderly as well as for those for whom the standard use of the bath or the shower poses a problem. Facilitates and assists during daily toilet activities Eg. at the sink, in the tub, showering.

Product information:

- » aluminium frame
- » seat made of easily washable polyethylene
- » U-cut seat to help wash the intimate parts
- » non-slip rubber tips
- » stable and robust design

Technical specification:

- » Seat dimensions: 41 x 48 cm
- » Min. height: 35 cm
- » Max. height: 45 cm
- » Weight: 2.1 kg
- » Weight capacity: 110 kg

Available colours: White





AR-204 Folding shower stool STEEL

General information:

Shower Stool as an assistive device is designed for the disabled with musculoskeletal system disorder, for elderly, as well as those for which using a standard bath or shower is problematic. Helps and assists during daily toilet activities eg. using washbasin, tub or shower stall.

Product information:

- » powder coated, steel frame
- » foldable frame (easy storage)
- » seat is made of an easily washable polyethylene
- » equipped with handles and drain holes
- » non-slip rubber tips
- » stable and robust design

Technical specs:

- » Seat dimensions: 50 x 30 cm
- » Height: 52 cm
- » Weight: 3.2 kg
- » Maximum load: 110 kg

Available colour: White



ARmedical®
HomeCare Line



ARmedical®
HomeCare Line

AR-205 Folding Shower Stool ALUMINIUM

General information:

Shower Stool as an assistive device is designed for the disabled with musculoskeletal system disorder, for elderly, as well as those for which using a standard bath or shower is problematic. Helps and assists during daily toilet activities eg. using washbasin, tub or shower stall.

Product information:

- » anodized aluminium frame
- » foldable frame (easy storage)
- » seat is made of an easily washable polyethylene
- » equipped with handles and drain holes
- » non-slip rubber tips
- » stable and robust design

Technical specs:

- » Seat dimensions: 50 x 30 cm
- » Min. height: 47 cm
- » Max. height: 50.5 cm
- » Weight: 2.4 kg
- » Maximum load: 110 kg

Available colour: White



AR-203 Shower chair with backrest

General information:

Shower chair as an assistive device is designed for the disabled with musculoskeletal system disorder, for elderly, as well as those for which using a standard bath or shower is problematic. Helps and assists during daily toilet activities eg. using washbasin, tub or shower stall.

Product information:

- » anodized aluminium frame
- » seat is made of an easily washable polyethylene
- » equipped with handles and drain holes
- » non-slip rubber tips
- » stable and durable construction with adjustable height

Technical specs:

- » Seat dimensions: 50 x 30 cm
- » Min. height: 36 cm
- » Max. height: 47 cm
- » Weight: 3 kg
- » Maximum load: 110 kg

Available colour: White



ARmedical®
HomeCare Line

**AR-208** Shower chair with U-cut**General information:**

Shower chair as a medical aid is designed for the disabled with locomotor disorders, the elderly as well as for those for whom the standard use of the bath or the shower poses a problem. Facilitates and assists during daily toilet activities eg. at the sink, in the tub, showering.

Product information:

- » anodized aluminium frame
- » seat and backrest made of easily washable polyethylene
- » U-cut seat to help wash the intimate parts
- » seat equipped with handles and drain holes
- » non-slip rubber tips
- » stable and robust design with height adjustment

Technical specification:

- » Seat dimensions: 48 x 33 cm
- » Min. height: 36 cm
- » Max. height: 56 cm
- » Weight: 2.7 kg
- » Weight capacity: 110 kg

Available colours: White**ARmedical®**
HomeCare Line**AR-207** Bathtub seat**General information:**




Bathtub bench as a medical designed for the disabled with locomotor disorders, the elderly, as well as people for whom the standard use of the bath poses a problem (sitting down and getting out). Facilitates and assists during daily toilet activities.

Product information:

- » aluminium frame
- » seat made of easily washable polyethylene
- » seat equipped with drain holes
- » adjustable handles width
- » stable and robust design

Technical specification:

- » Seat dimensions: 41 x 22 cm
- » Recommended width of the bathtub: 52 – 61 cm
- » Weight: 1.5 kg
- » Weight capacity: 110 kg

**ARmedical®**
HomeCare Line**ARmedical®**
HomeCare Line**NEW soft, ergonomic
non-slip handle****AR-206** Bathtub bench with handle**General information:**

Bathtub bench as an assistive device is designed for the disabled with musculoskeletal system disorder, for elderly, as well as those for which using a standard tub is problematic (sitting down and standing up). Facilitates and assists during daily toilet routine..

Product information:

- » seat made of an easily washable polyethylene
- » seat equipped with drain holes and place for soap
- » handle to help while standing up
- » adjustable expanders
- » stable and robust design

Technical specs:

- » Seat dimensions: 70 x 30 cm
- » Seat height: 5 cm
- » For the tub width: 40 – 62 cm
- » Weight: 2.7 kg
- » Maximum load: 160 kg

Available handle colours: Blue or Green

Medical products

Hot & Cold Packs



Reusable cold-hot compress with wide range of applications. As thermal therapy causes dilation of blood vessels at the site of action, raising the temperature of tissues, increasing blood flow and stimulation of cellular metabolism, modulates the treatment of inflammatory conditions. As cold therapy it is soothing in case of headache, toothaches. Promotes the inhibition of edema in case of injuries (contusions, bruises, sprains, etc.).

- Perfect as a soothing headaches and contusions
- Ideal for children and adults as an analgesic and anti-edema aid

Application:

Warm Compress

application as a warm compress in the case of:

- » rheumatic pain
- » menstrual pain
- » chronic pain (eg. lumbago),
- » painful muscle spasms

Cold Compress


application as a cold compress in the case of:

- » injuries (contusions, sprains, etc.)
- » headaches
- » toothaches
- » swelling



CHP-1010

Cold-Hot Compress
"Mini", 10x10cm



CHP-1025

Cold-Hot Compress
"Classic", 10x25cm



CHP-2030

Cold-Hot Compress
"Maxi", 20x30cm

ARmedical®
Hot & Cold Therapy

A series of horizontal dotted lines spanning the width of the page, providing a template for text entry.



s.1	MFP-5030	Profled orthopaedic pillow "EXCLUSIVE DREAM"
s.2	MFP-5030MV	Profled orthopaedic pillow with massage & ventilation "EXCLUSIVE DREAM"
s.3	MFP-5535 MFP-5030BF	Profled orthopaedic pillow "EXCLUSIVE DREAM" Profled orthopaedic pillow "BAMBOO DREAM"
s.4	MFP-5030BV MFP-5232	Profled orthopaedic pillow with massage & ventilation "Bamboo Dream" Profled orthopaedic pillow "BAMBOO DREAM"
s.5	MFP-5030AF MFP-5535AF	Profled orthopaedic pillow "AloeVera Dream" Profled orthopaedic pillow "AloeVera Dream"
s.6	MFP-5433EP	Double-sided profled orthopaedic pillow "Ergonomic Dream"
s.7	MFP-6035	Profled orthopaedic pillow "Silent Dream"
s.8	TFP-5030TP TFP-6040	Traditional orthopaedic pillow "EXCLUSIVE DREAM" Traditional orthopaedic pillow "EXCLUSIVE DREAM"
s.9	MFP-4026	Profled orthopaedic pillow "BABY DREAM"
s.10	MTP-2828 MTP-3230	Orthopaedic pillow "EXCLUSIVE TRAVEL" Orthopaedic pillow "EXCLUSIVE TRAVEL"
s.11	MFP-3433 MFP-5344	Orthopaedic lumbar pillow "EXCLUSIVE SUPPORT" Orthopaedic lumbar pillow "EXCLUSIVE SUPPORT"
s.12	MFP-2520 MFP-4220	Orthopaedic knee "separator" pillow "EXCLUSIVE SUPPORT" Orthopaedic lumbar pillow „EXCLUSIVE SUPPORT"
s.13	DSC-33	Dynamic seat cushion "DYNAPAD"
s.14	DSC-35 DSC-37	Dynamic seat cushion "DYNAPAD MAX" Dynamic seat cushion "DYNAPAD BALANCE"
s.15	GM-55 to 85 ABS-55 to 85	Rehabilitation ball "CLASSIC" ABS Rehabilitation ball
s.16	BMB-55 to 15 MB-5 do 10 TMB-20	Massage rehabilitation ball Squeeze rehabilitation ball with spikes Massage rehabilitation ball
s.17	RLB PSB-50	Rehabilitation ball MIDI REH ABS Rehabilitation ball "PEANUT"
s.18	HB2-45 to 55 HB1-45 to 55	Rehabilitation ball "HOPPER" with horn handles Rehabilitation ball "HOPPER" with handle
s.19	AR-920 AR-950	Anti Decubitis Air Bubble Mattress & Pump System Anti Decubitis Air Tubing Mattress & Pump System
s.20	AR-018 AR-019	Pedal exerciser Pedal exerciser – folding with meter
s.21	AR-010 AR-013	Elbow crutch – aluminium Aluminium T-shape cane
s.22	AR-014 AR-015	Aluminium T-shape cane with soft handgrip AR-015 Aluminium T-shape cane – foldable
s.23	AR-016 AR-017	Aluminium T-shape Tripod Cane Aluminium T-shape Quad Cushion Cane with soft handgrip
s.24	AR-017W AR-001	Aluminium T-shape High Quad Cushion Cane Aluminium folding walker. Rolling – stable. Single H-shaped frame
s.25	AR-002 AR-003	Aluminium folding walker Rolling – stable. Double H-shaped frame Aluminium folding walker with two 5-inch wheels
s.26	AR-008 AR-009	Aluminium folding walker – 3 functional Aluminium folding walker with ball – Two-wheels
s.27	AR-020 AR-021	Two-wheel rehabilitation walker STANDARD Two-wheel walker DELUXE
s.28	AR-004 AR-005	3-wheel Rollator with storage bag and basket 4-wheel steel rehab Rollator CLASSIC
s.29	AR-006 AR-007	4-wheel aluminium rehab Rollator PRESTIGE 4-wheel aluminium rehab Rollator DYNAMIC
s.30	AR-022	Four-wheel aluminium rollator with footrest COMFORT
s.31	AR-300	Aluminium wheelchair ERGONOMIC
s.32	AR-330	Aluminium wheelchair DYNAMIC
s.33	AR-350	Aluminium wheelchair PRESTIGE
s.34	AR-101 AR-110/AR-115	Folding – steel toilet chair Raised toilet seat with lid
s.35	AR-201 AR-202 AR-209	Rounded shower stool Shower chair without backrest Shower stool with U-cut
s.36	AR-204 AR-205 AR-203	Folding shower stool STEEL Folding Shower Stool ALUMINIUM Shower chair with backrest
s.37	AR-208 AR-207 AR-206	Shower chair with U-cut Bathtub seat Bathtub bench with handle
s.38	CHP-1010 CHP-1025 CHP-2030	Cold-Hot Compress "Mini", 10x10cm Cold-Hot Compress "Classic", 10x25cm Cold-Hot Compress "Maxi", 20x30cm

More steadily through life!



Please feel invited to cooperate with us!

ARmedical®

41-303 Dąbrowa Górnicza, ul. Rozdzieńskiego 19
tel. fax: +48 32 261 71 36; fax: +48 32 268 63 05
e-mail: biuro@armedical.pl
www.armedical.pl
www.rehabilitacyjnepilki.pl
www.ortopedycznepoduszki.pl
www.lekkigips.pl

Regional Distributor

ARmedical®

is registered trademark of ARMEDICAL S.C.

CE

ARmedical® products are reported to the Office for Registration of Medicinal Products as a class I medical device in accordance with the requirements of the Medical Devices Act of 20.05.2010r. and Medical Directive (93/42/EEC) concerning medical devices. Manufacturer of medical devices ARmedical S.C. issued the appropriate EC Declaration of Conformity.

This Product Catalogue cannot be considered as a trade offer within the meaning of the Civil Code (Art. 71kc) is for informational purposes only.

All information contained in this catalogue for the products manufactured at the request, imported and distributed by the ARmedical company that can be entered into the system of reimbursement from public funds are for informational purposes only and are not an endorsement in accordance within the meaning of Art. 49.ust.3 Act on medicine reimbursement, special purpose food agents and medical products dated 12.05.2011.